

## REGULAR LAB 07d: Limiting Reactant & % Yield

**Aim** To investigate the stoichiometry of the reaction between sodium hydrogen carbonate and ethanoic acid. The reaction is;



**Apparatus** Four 125 mL Erlenmeyer flasks, graduated cylinders, filter papers, electronic balance

**Chemicals** Sodium hydrogen carbonate, ethanoic acid

### **Method**

1. Weigh five, 3.50 g samples of sodium hydrogen carbonate, onto five different filter papers.
2. Using a graduated cylinder, measure 10.0 mL of ethanoic acid. Pour the acid into a clean, 125 mL Erlenmeyer flask. Rinse the graduated cylinder with a **small** amount of water and pour this rinse water into the flask. Measure the mass of the flask and contents.
3. Pour one of the 3.50 g samples of sodium hydrogen carbonate into the flask and swirl to mix the contents **until no more bubbling occurs** (this may take a few minutes). **Do not allow any of the contents to splash out.**
4. **When all bubbling has ended**, determine the mass of the flask and contents and record in the results table. Pour out the contents, wash and dry the flask.
5. Repeat steps #2 through #4, four more times, using 30.0, 50.0, 70.0 and 90.0 mL of ethanoic acid respectively instead of 10.0 mL.

## Results

		EXPERIMENT				
		1	2	3	4	5
<b>A</b>	Mass of NaHCO <sub>3</sub> in g	3.50	3.50	3.50	3.50	3.50
<b>B</b>	Molar mass of NaHCO <sub>3</sub> in g mol <sup>-1</sup>	84.0	84.0	84.0	84.0	84.0
<b>C</b>	Moles of NaHCO <sub>3</sub> in mols					
<b>D</b>	Volume of CH <sub>3</sub> COOH in mL	10.0	30.0	50.0	70.0	90.0
<b>E</b>	Moles* of CH <sub>3</sub> COOH in mols					
<b>F</b>	Excess Reagent					
<b>G</b>	Limiting Reagent					
<b>H</b>	Theoretical Mass Loss in g					
<b>I</b>	Mass of flask + acid + water in g					
<b>J</b>	Mass of flask + acid + rinse water + NaHCO <sub>3</sub> in g					
<b>K</b>	Final mass of flask + contents in g					
<b>L</b>	Actual Mass Loss in g					
<b>M</b>	% Yield					

\* The ethanoic acid is a solution that has a concentration of 0.837M (or 0.837 moles per liter of solution). Moles of the acid can be calculated using moles = (concentration) x (volume in L).

**Conclusion/Calculations:**

1. Make a graph of your data, using mass loss on the y-axis with moles of ethanoic acid on the x-axis. Use (0,0) as the origin.
2. Why did the flasks lose mass?
3. Is there a point at which adding more ethanoic does not increase the mass loss? Explain.



4. Is there a limit to the quantity of gas that can be produced from 3.50 grams of sodium hydrogen carbonate? Explain.